


available at [www.sciencedirect.com](http://www.sciencedirect.com)journal homepage: [www.elsevier.com/locate/peptides](http://www.elsevier.com/locate/peptides)

# Identification of PVK/CAP2b neuropeptides from single neurohemal organs of the stable fly and horn fly via MALDI-TOF/TOF tandem mass spectrometry

Ronald J. Nachman<sup>a,\*</sup>, William K. Russell<sup>b</sup>, Geoffrey M. Coast<sup>c</sup>, David H. Russell<sup>b</sup>, J. Allen Miller<sup>d</sup>, Reinhard Predel<sup>a,e</sup>

<sup>a</sup> Areawide Pest Management Research, Southern Plains Agricultural Research Center, USDA, 2881 F/B Road, College Station, TX 77845, USA

<sup>b</sup> The Laboratory for Biological Mass Spectrometry, Department of Chemistry, Texas A&M University, College Station, TX 77843, USA

<sup>c</sup> School of Biological and Chemical Sciences, Birkbeck, University of London, London WC1E 7HX, UK

<sup>d</sup> Knippling-Bushland U.S. Livestock Entomology Research Center, USDA, Kerrville, TX 78028, USA

<sup>e</sup> Saxon Academy of Sciences, Research Group Jena, Erbertstrasse 1, Jena 07743, Germany

## ARTICLE INFO

### Article history:

Received 14 July 2005

Accepted 28 July 2005

Published on line 2 December 2005

### Keywords:

MALDI-TOF/TOF mass spectrometry

Insect

Neuropeptide

Periviscerokinin

CAP2b

*Stomoxys calcitrans*

*Haematobia irritans*

Peptidomics

## ABSTRACT

MALDI-TOF/TOF tandem mass spectrometry has been applied to determine the complete sequences of the PVK/CAP2b neuropeptides in the stable fly *Stomoxys calcitrans* and horn fly *Haematobia irritans*, insect pests of livestock. This peptidomic analysis of single neurohemal organ preparations allows the unambiguous assignment of internal Leu/Ile positions not distinguishable by previous mass spectrometric techniques. The sequences are as follows: Stoca-PVK/CAP2b-1, AGGASGLYAFPRVa; Stoca-PVK/CAP2b-2, NAKLYPVPRVa; and Haeir-PVK/CAP2b-1, AGGASGLYAFPRVa; Haeir-PVK/CAP2b-1, NAKLYPMPRVa. Both Stoca-PVK/CAP2b-1 and -2 stimulate Malpighian tubule fluid secretion in the stable fly, with EC<sub>50</sub> values between 3 and 11 nM. The identification of these novel neuropeptides adds to our knowledge of the peptidomes of flies, and can aid in the development of neuropeptide-based control strategies of these insect pests.

Published by Elsevier Inc.

## 1. Introduction

Neuropeptides are important messenger molecules that occur in a great variety of forms and are implicated in the regulation of critical physiological processes such as diuresis, digestion, development and reproduction [6,7] in insects. In the past several years, new developments in matrix-assisted laser desorption-time-of-flight mass spectrometry (MALDI-TOF MS) have afforded very sensitive de novo sequencing of peptides

via direct analysis of single neurosecretory organs or nerves, including those of insects [13,18,19] via the post-source decay (PSD) technique. Alternatively, electrospray ionization (ESI) coupled with tandem MS has resulted in the identification of a number of novel neuropeptides [2,20,22,24,25]. The amount of material needed for ESI-MS experiments is usually larger than the amount necessary for MALDI-MS, since the peptides have to be extracted from the tissue prior to analysis. Both techniques, however, alleviate the need for large numbers

\* Corresponding author. Tel.: +1 979 260 9315; fax: +1 979 260 9377.

E-mail address: [nachman@tamu.edu](mailto:nachman@tamu.edu) (R.J. Nachman).

0196-9781/\$ – see front matter. Published by Elsevier Inc.

doi:10.1016/j.peptides.2005.07.022

Report Documentation Page				Form Approved OMB No. 0704-0188	
Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.					
1. REPORT DATE <b>JUL 2005</b>		2. REPORT TYPE		3. DATES COVERED <b>00-00-2005 to 00-00-2005</b>	
4. TITLE AND SUBTITLE <b>Identification of PVK/CAP2b neuropeptides from single neurohemal organs of the stable fly and horn fly via MALDI-TOF/TOF tandem mass spectrometry</b>				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) <b>U.S. Department of Agriculture,Areawide Pest Management Research,2881 F/B Road,College Station,TX,77845</b>				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT <b>Approved for public release; distribution unlimited</b>					
13. SUPPLEMENTARY NOTES					
14. ABSTRACT <b>see report</b>					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT <b>Same as Report (SAR)</b>	18. NUMBER OF PAGES <b>6</b>	19a. NAME OF RESPONSIBLE PERSON
a. REPORT <b>unclassified</b>	b. ABSTRACT <b>unclassified</b>	c. THIS PAGE <b>unclassified</b>			

of specimens and the time-consuming and expensive efforts required to isolate and determine the primary structure of neuropeptides via traditional chromatographic and chemical sequencing techniques [4,8,28]. These MS techniques have also made it possible to make detailed comparisons of the peptide patterns or profiles (i.e., the peptidomes) of closely related insect species [1,3,24]. MS analysis of insect neuropeptides has failed in the past to distinguish between the isomers Leu and Ile, which have identical masses. These earlier studies were limited to low energy fragmentations of the ion of interest. A primary limitation of PSD peptide sequencing is that the internal energies of the  $[M + H]^+$  ions are not sufficient to yield the side chain cleavages necessary to distinguish Leu and Ile. However, recent innovations in MALDI-TOF MS have allowed analysis of high-energy collision-induced dissociation of the parent ions of peptides that reveal unique mass patterns for the sidechains of Leu and Ile [11,12]. Indeed, in a recent study we have demonstrated the utility of MALDI-TOF/TOF tandem MS in distinguishing between Leu and Ile in neuropeptides from single neurohemal preparations of insects, specifically PVK/CAP2b (periviscerokinin/cardioacceleratory peptide 2b) sequences from the housefly (*Musca domestica*) and flesh fly (*Neobellieria bullata*) [15].

In this study, we utilize MALDI-TOF/TOF tandem MS to undertake the first identification of the sequences of PVK/CAP2b neuropeptides, including the unambiguous assignment of Leu/Ile positions, from single neurohemal organ preparations of adults of the stable fly *Stomoxys calcitrans* and horn fly *Haematobia irritans*, important livestock pests. PVK/CAP2bs are typical of the abdominal neurohemal system of insects [29], usually stored in abdominal perisymphatic organs. The PVK/CAP2b class of neuropeptides has been shown to elicit both myotropic activity and stimulation of Malpighian tubule fluid secretion in insects (see [29]), physiological processes critical to survival. We report on the diuretic activity of the two PVK/CAP2b sequences native to the stable fly in a stable fly Malpighian tubule fluid secretion bioassay.

## 2. Materials and methods

### 2.1. Insects

Stable flies (*S. calcitrans*) and horn flies (*H. irritans*) were obtained as pupae from the Knipling-Bushland U.S. Livestock Insects Laboratory in Kerrville, TX. Pupae were kept in cages at 26 °C. Newly emerged stable and horn flies were fed a beef blood meal twice daily for two weeks and provided with sucrose and water ad libitum. Male and female flies <10 days post-emergence were used in all experiments. Stable flies used in fluid secretion assays were obtained as pupae from Professor Mike Lehane (Liverpool School of Tropical Medicine). Adult flies were held at 27 °C and fed a 5% solution of sucrose.

### 2.2. Matrix-assisted laser desorption ionization tandem time-of-flight mass spectrometry (MALDI-TOF/TOF MS)

Dissection and sample preparation were performed as previously described [14,18]. MALDI analysis was performed

on the ABI 4700 proteomics analyzer (Applied Biosystems, Framingham, MA) [12]. Due to the nature of the samples all acquisitions were taken in manual mode. Initially the instrument was operated in reflectron mode, in order to determine the parent masses. The laser intensity was set just above the threshold required to ionize the neuropeptides. For the tandem MS experiments, the acceleration voltage applied was 1 kV in all cases, and the laser intensity was increased by 10%. The number of laser shots used to obtain a spectrum varied from 500 to 5000, depending on signal quality. In order to change the net amount of collisions to the primary ions in the collision induced dissociation (CID) experiment, the collision cell gas (atmospheric air) pressure was increased. All three gas pressures settings available (none, medium and high) were employed. The instrument was operated in post-source decay mode when no collision gas is used. The true pressure within the collision cell cannot be measured. The fragmentation patterns from these three different settings were used to determine the sequence of the peptide. The fragmentation data obtained in these experiments was handled using the Applied Biosystems Data Explorer<sup>®</sup> software package.

### 2.3. Peptide synthesis

Stoca-PVK-1/Haair-PVK-1 (AGGASGLYAFPRVa), Stoca-PVK-2 (NAKLYPVPRVa), and Haair-PVK-2 (NAKLYPMPRVa) were synthesized via Fmoc methodology on Rink Amide resin (Novabiochem, San Diego, CA) using Fmoc protected amino acids (Advanced Chemtech, Louisville, KY) on an ABI 433A peptide synthesizer (Applied Biosystems, Foster City, CA) under previously described conditions [14]. Crude products were purified on a Waters C<sub>18</sub> Sep Pak cartridge and a Delta Pak C<sub>18</sub> reverse-phase column (8 mm × 100 mm, 15 μm particle size, 100 Å pore size) on a Waters 510 HPLC controlled with a Millennium 2010 chromatography manager system (Waters, Milford, MA) with detection at 214 nm at ambient temperature. Solvent A = 0.1% aqueous trifluoroacetic acid (TFA); Solvent B = 80% aqueous acetonitrile containing 0.1% TFA. Conditions: Initial solvent consisting of 20% B was followed by the Waters linear program to 100% B over 40 min; flow rate, 2 ml/min. Delta-Pak C-18 retention times: Stoca-PVK-1/Haair-PVK-1 (AGGASGLYAFPRVa), 9.0 min; Stoca-PVK-2 (NAKLYPVPRVa), 6.25 min; and Haair-PVK-2 (NAKLYPMPRVa), 10.95 min. The peptides were further purified on a Waters Protein Pak I125 column (7.8 mm × 300 mm) (Milligen Corp., Milford, MA). Conditions: Flow rate: 2.0 ml/min; isocratic with solvent = 80% acetonitrile made to 0.01% TFA. WatPro retention times: Stoca-PVK-1/Haair-PVK-1 (AGGASGLYAFPRVa), 9.25 min; Stoca-PVK-2 (NAKLYPVPRVa), 8.75 min; and Haair-PVK-2 (NAKLYPMPRVa), 9.0 min. Amino acid analysis was carried out under previously reported conditions [14] and used to quantify the peptide and to confirm identity, leading to the following analyses: Stoca-PVK-1/Haair-PVK-1 (AGGASGLYAFPRVa): A[2.8], F[1.0], G[2.7], L[1.0], P[1.1], R[1.0], S[0.9], V[1.1], Y[1.0]; Stoca-PVK-2 (NAKLYPVPRVa), A[1.0], K[0.9], L[1.0], N[0.9], P[1.9], R[1.0], V[1.9], Y[1.0]; and Haair-PVK-2 (NAKLYPMPRVa), A[1.0], K[1.0], L[1.0], M[1.0], N[1.0], P[2.0], R[1.0], V[1.0], Y[0.9]. The identity of the peptide analogs was confirmed via MALDI-TOF-MS on a Kratos Kompact Probe MALDI-TOF MS machine (Kratos Analytical, Ltd., Manchester,

UK) with the presence of the following molecular ions ( $M + H^+$ ): Stoca-PVK-1/Haeir-PVK-1 (AGGASGLYAFPRVa), 1266.3 [ $M + H^+$ ]; Stoca-PVK-2 (NAKL $\underline{Y}$ VPVRVa), 1156.6 [ $M + H^+$ ]; and Haeir-PVK-2 (NAKL $\underline{Y}$ PMPRVa), 1188.9 [ $M + H^+$ ].


## 2.4. Isolated stable fly Malpighian tubule preparations

Fluid secretion from isolated stable fly Malpighian tubules was measured using the “Ramsay assay” as previously described for housefly tubules [8]. Tubules were removed from 4- to 7-day-old adult flies of both sexes. Flies were dissected under Musca saline [8] and both anterior and posterior tubules were transferred to small (10  $\mu$ l) drops of bathing fluid (a 1:1 mixture of saline and Schneider’s *Drosophila* medium) beneath water-saturated liquid paraffin in a Sylgard<sup>TM</sup> lined Petri dish. The tubules were allowed to equilibrate for 1 h before being challenged with test peptides. Urine escaped from the cut end of the tubule, which was pulled out into the liquid paraffin. Drops of urine were collected at 15 min intervals and their diameter ( $d$ ) measured as they rested on the Sylgard base of the Petri dish using a Wild digital (MMS235) eyepiece micrometer. Urine volume was calculated as  $\pi d^3/6$  and the rate of secretion obtained by dividing the secreted volume by the collection period. Data were normalized by expressing the increase in fluid secretion as a percentage of the response to 10 nM Musdo-K (also native to the stable fly [13]), which was added to all tubules at the end of each experiment. Dose-response curves were prepared using the computer program GraphPad Prism version 4.02 (GraphPad Software, San Diego, CA).


## 3. Results

### 3.1. Determination of PVK/CAP2b sequences in abdominal dorsal sheaths of flies via MALDI-TOF/TOF mass spectrometry

Direct analysis of abdominal dorsal sheath tissues from the stable fly *S. calcitrans* and horn fly *H. irritans* were conducted via MALDI-TOF/TOF MS. Illustrated in Fig. 1 are the initial MALDI-TOF spectra of the preparations of these two flies that feature the parent ions of the PVK/CAP2b peptides. High energy CID of the PVK/CAP2b peptides reveal unique patterns for the sidechains of Leu and Ile [11,12,15]. As illustrated in Fig. 2, fragments of native Stoca-PVK-2 of the stable fly include a prominent ‘w7a’ fragment ion at mass 788.3 Da (Fig. 2), indicative of Leu in this peptide. Indeed, the spectrum of the synthetic version of Stoca-PVK-2 containing Leu at position 7 taken under the same conditions was essentially identical. In contrast, if Stoca-PVK-2 contained Ile, the mass spectrum under conditions of high gas would reveal a different mass for the ‘w7a’ fragment, along with two diagnostic satellite ions, a ‘v-ion’ and a ‘wb-ion’. Thus, Stoca-PVK-2 can be unambiguously assigned the sequence NAKLYPVPRVa (Table 1). Using the same single organ preparations, the sequences of Stoca-PVK-1, Haeir-PVK-1, and Haeir-PVK-2 could also be unambiguously determined under conditions of high gas pressure to be AGGTSGLYAFPRVa, AGGTSGLYAFPRVa, and NAKLYPMPRVa (Table 1), respectively.



**Fig. 1** – MALDI-TOF mass spectra of preparations of the posterior abdominal dorsal sheath of the ventral nerve cord of the stable fly *S. calcitrans* (A) and the horn fly *H. irritans* (B). The dorsal sheath represents a neurohemal release site homologous to perisymphatic organs which become incorporated into the dorsal ganglionic sheath during the metamorphosis of cycloraphous flies. The  $MH^+$  values are as follows: for Stoca-PVK/CAP2b-1 (1264.70), Stoca-PVK/CAP2b-2 (1155.72), Haeir-PVK/CAP2b-1 (1264.68), and Haeir-PVK/CAP2b-2 (1187.67).



**Fig. 2** – A MALDI-TOF/TOF tandem mass spectrum of native Stoca-PVK-2 under conditions of high gas pressure (arrow identifies the ‘w7a’ fragment ion at 788.3). The mass spectra of the native and synthetic Leu variant are essentially identical.

**Table 1 – Amino acid sequences of PVK/CAP2b peptides native to the stable fly (*S. calcitrans*) and horn fly (*H. irritans*) determined by MALDI-TOF-TOF tandem mass spectrometry compared with previously determined sequences from other flies**

Species	PVK/CAP2b-1	PVK/CAP2b-2
<i>S. calcitrans</i>	AGGASGLYAFPRVa	NAKLYPVPRVa
<i>H. irritans</i>	AGGASGLYAFPRVa	NAKLYPMRPRVa
<i>M. domestica</i>	AGGTSGLYAFPRVa	ASLFNAPRVa [15]
<i>N. bullata</i>	NGGTSGLYAFPRVa	AGLIVYPR[L/I]a <sup>a</sup> [15]
<i>D. melanogaster</i>	GANMGLYAFPRVa	ASGLVAFPRVa [27]
<i>An. gambiae</i>	GPTVGLYAFPRVa	QGLVPFPRVa [17]


<sup>a</sup> MALDI-TOF/TOF tandem MS cannot distinguish between Leu and Ile at a C-terminal position [14].

### 3.2. Fluid secretion activity of PVK/CAP2b sequences on stable fly Malpighian tubules

Synthetic replicates of the two PVK/CAP2b sequences native to the stable fly were evaluated for diuretic activity on isolated Malpighian tubules. As can be seen in Fig. 3 and Table 2, Stoca-PVK-1 and Stoca-PVK-2 stimulated stable fly Malpighian tubule fluid secretion with observed EC<sub>50</sub> values of 3.4 and 11.0 nM, respectively (Table 2). The peak response to both peptides was ~25% of that obtained after the addition of 10 nM Musdo-K (also native to the stable fly [13]), and represented an 80% increase over the basal rate of fluid secretion.

## 4. Discussion

PVK/CAP2b peptides have been found in the stable fly *S. calcitrans* and horn fly *H. irritans* via MALDI-TOF MS analysis of preparations of neurohemal tissues [24]. PVK/CAP2b peptides of insects are typical of neurosecretory neurons in the



**Fig. 3 – Dose–response curves for Stoca-PVK/CAP2b-1 (solid line and symbols) and -2 (dashed line and open symbols). Data are expressed as the mean  $\pm$  S.E.M for six replicates.**

**Table 2 – Fluid secretion activity of native Stomoxys-PVK/CAP2b-1 and -2 on stable fly (*S. calcitrans*) Malpighian tubules**

Peptide	Sequence	Fluid secretion <sup>a</sup> [EC <sub>50</sub> ] (nM)
Stoca-PVK/CAP2b-1	AGGASGLYAFPRVa	3.4 (0.6–18.9)
Stoca-PVK/CAP2b-2	NAKLYPVPRVa	11.0 (3.8–31.6)

<sup>a</sup> Values in parentheses, represent the 95% confidence limit (CL).

abdominal ganglia and are accumulated in perisymphatic organs until release. Larval perisymphatic organs of cyclorhous Diptera, however, become incorporated into the dorsal ganglionic sheath [16] during the metamorphosis. The abdominal dorsal sheath, which was dissected in this study, therefore contains relatively large amounts of peptidergic neurohormones. The molecular ions and incomplete sequences observed in the MALDI-TOF/TOF mass spectrometric studies operated in reflectron mode (Fig. 1) on the stable fly and horn fly were, Stoca-PVK-1 (AGGASG[L/I]YAFPRVa; *m/z* 1264.70) and Stoca-PVK-2 (NAK[L/I]YPVPRVa; *m/z* 1155.72); and Haeir-PVK-1 (AGGASG[L/I]YAFPRVa; *m/z* 1264.68) and Haeir-PVK-2 (NAK[L/I]YPMRPRVa; *m/z* 1187.67). When the peptides were fragmented under conditions of high collision energy, the collision-induced fragments reveal distinct side-chain fragmentations. Stoca-PVK-2 demonstrated a fragmentation pattern indicative of a Leu in position 7 (Fig. 2), and this same fragmentation pattern was observed in mass spectra of a synthetic replicate containing Leu. Thus, Stoca-PVK-1 can be unambiguously assigned the sequence NAKLYPVPRVa (Table 1). Using the single organ preparations, the sequences of Stoca-PVK-1, Haeir-PVK-1 and Haeir-PVK-2 could also be unambiguously determined under conditions of high collision energy to be AGGASGLYAFPRVa, AGGASGLYAFPRVa and NAKLYPMRPRVa (Table 1), respectively. Thus, the PVK/CAP2b-1 sequences of the stable fly and horn fly are identical, and the PVK/CAP2b-2 sequences differ only by a single residue in the fourth position from the C-terminus (V versus M).

It is clear that the Leu at the position located seven residues from the C-terminus in these four PVK/CAP2b sequences from stable fly and horn fly is conserved within and across species (Table 1). Leu at this specific position is typical of other PVK/CAP2b of insects that were sequenced by Edman degradation in earlier studies [5,9,20,21,23] or for which genes have been published already [10,26,27]. In general, PVK/CAP2b sequences are well conserved throughout the six species of flies that have been studied to date (Table 1).

Both the stable fly peptides stimulated Malpighian tubule fluid secretion, with Stoca-PVK/CAP2b-1 being slightly more potent than Stoca-PVK/CAP2b-2, although the difference was not significant. The maximum response was equivalent to an 80% increase in the basal rate of secretion, which is consistent with data for other dipteran insects [17], and is about 25% of the response of stable fly tubules to 10 nM Musdo-K.

In summary, MALDI-TOF/TOF tandem MS has been used to identify the sequences for the PVK/CAP2b neuropeptides in the stable fly and horn fly via direct analysis of nerve tissue. This analysis includes an unambiguous assignment of the monoisotopic residues Leu versus Ile. The work adds to our knowledge of the peptidomes of flies; and the identification of



the specific structures of the PVK/CAP2b neuropeptides, implicated in the regulation of diuretic and myotropic processes, may aid in the development of mimetic analogs capable of disrupting these critical physiological processes in pest flies.

## Acknowledgements

This study was supported by a Binational Agricultural Research and Development Grant (BARD #IS-3356-02) (RJN), a grant from Deutsche Forschungsgemeinschaft (Predel 595/6-1), a Collaborative Research Grant (#LST.CLG.979226) from the North Atlantic Treaty Organization (NATO) (RJN, GMC), a grant from the USDA/DOD DWFP Initiative (#0500-32000-001-01R) (RJN), and grants provided to the Laboratory for Biological Mass Spectrometry at Texas A&M University, College Station: R01 RR019587 NIH and DE-FG02-04ER15520 DOE. In addition, we acknowledge the capable technical assistance of Allison Strey, Pawel Zubrzak, and Nan Pryor of the Areawide Pest Management Research Unit, Southern Plains Agricultural Research Center.

## REFERENCES

- [1] Audsley N, Weaver RJ. Identification of neuropeptides from brains of larval *Manduca sexta* and *Lacanobia oleracea* using MALDI-TOF mass spectrometry and post-source decay. *Peptides* 2003;24:1465–74.
- [2] Baggerman G, Cerstiaen A, De Loof A, Schoofs L. Peptidomics of the larval *Drosophila melanogaster* central nervous system. *J Biol Chem* 2002;277:40368–74.
- [3] Clynen E, Baggerman G, Huybrechts J, Vanden Bosch L, De Loof A, Schoofs L. Peptidomics of the locust corpora allata: identification of novel pyrokinins (FXPRamide). *Peptides* 2003;24:1493–500.
- [4] Clynen E, Baggerman G, Veelaert A, Cerstiaens A, Van der Horst D, Harthoorn L, et al. Peptidomics of the pars intercerebralis-corpora cardiaca complex of the migratory locust, *Locusta migratoria*. *Eur J Biochem* 2001;268:1929–39.
- [5] Clynen E, Huybrechts J, De Loof A, Schoofs L. Mass spectrometric analysis of the perisymphathetic organs in locusts: identification of novel periviscerokinins. *Biochem Biophys Res Commun* 2003;300:422–8.
- [6] Coast GM. The neuroendocrine regulation of salt and water balance in insects. *Zoology (ZACS)* 2001;103:179–88.
- [7] Gäde G. The explosion of structural information on insect neuropeptides. In: Herz W, Kirby GW, Moore RE, Steglich W, Tamm C, editors. *Progress in the chemistry of organic natural products*. Viena: Springer-Verlag; 1997. p. 1–128.
- [8] Holman GM, Nachman RJ, Coast GM. Isolation, characterization and biological activity of a diuretic myokinin neuropeptide from the housefly, *Musca domestica*. *Peptides* 1999;20:1–10.
- [9] Huesmann GR, Cheung CC, Loi PK, Lee TD, Swiederek KM, Tublitz NJ. Amino acid sequence of CAP<sub>2b</sub>, an insect cardioacceleratory peptide from the tobacco hawkmoth *Manduca sexta*. *FEBS Lett* 1995;371:311–4.
- [10] Loi PK, Tublitz NJ. Sequence and expression of the CAPA/CAP2b gene in the tobacco hawkmoth, *Manduca sexta*. *J Exp Biol* 2004;207:3681–91.
- [11] Macht M, Asperger A, Deininger SO. Comparison of laser-induced dissociation and high-energy collision-induced dissociation using matrix-assisted laser desorption/ionization tandem time-of-flight (MALDI-TOF/TOF) for protein and peptide identification. *Rapid Commun Mass Spectrom* 2004;18:2093–105.
- [12] Medzihradsky KF, Campbell JM, Baldwin MA, Falick AM, Juhasz P, Vestal ML, et al. The characteristics of peptide collision-induced dissociation using a high-performance MALDI-TOF/TOF tandem mass spectrometer. *Anal Chem* 2000;72:552–8.
- [13] Nachman RJ, Coast GM, Tichy SE, Russell DH, Miller JA, Predel R. Occurrence of insect kinins in the flesh fly, stable fly and horn fly—mass spectrometric identification from single nerves and diuretic activity. *Peptides* 2002;23:1885–94.
- [14] Nachman RJ, Isaac RE, Coast GM, Holman GM. Aib-containing analogues of the insect kinin neuropeptide family demonstrate resistance to an insect angiotensin-converting enzyme and potent diuretic activity. *Peptides* 1997;18:53–7.
- [15] Nachman RJ, Russell DH, Coast GM, Russell DH, Predel R. Mass spectrometric assignment of Leu/Ile in neuropeptides from single neurohemal organ preparations of insects. *Peptides* 2005;26:2151–6.
- [16] Nässel DR, Ohlsson LG, Cantera R. Metamorphosis of identified neurons innervating thoracic neurohemal organs in the blowfly—transformation of cholecystokinin like immunoreactive neurons. *J Comp Neurol* 1988;267:343–56.
- [17] Pollock VP, McGettigan J, Cabrero P, Maudlin IM, Dow JAT, Davies SA. Conservation of capa peptide-induced nitric oxide signaling in Diptera. *J Exp Biol* 2004;207:4135–45.
- [18] Predel R, Gäde G. Identification of the abundant neuropeptide from abdominal perisymphathetic organs of locusts. *Peptides* 2002;23:621–8.
- [19] Predel R, Gäde G. Peptidomics of neurohemal organs from species of the cockroach family Blattidae: how do neuropeptides of closely related species differ? *Peptides* 2005;26:3–9.
- [20] Predel R, Kellner R, Baggerman G, Steinmetzer T, Schoofs L. Identification of novel periviscerokinins from single neurohaemal release sites in insects—MS/MS fragmentation complemented by Edman degradation. *Eur J Biochem* 2000;267:3869–74.
- [21] Predel R, Linde D, Rapus J, Vettermann S, Penzlin H. Periviscerokinin (Pea-PVK): a novel myotropic neuropeptide from the perisymphathetic organs of the American cockroach. *Peptides* 1995;16:61–6.
- [22] Predel R, Neupert S, Wicher D, Gundel M, Roth S, Derst C. Unique accumulation of neuropeptides in an insect: FMRamide related peptides in the cockroach, *Periplaneta americana*. *Eur J Neurosci* 2004;20:1499–513.
- [23] Predel R, Rapus J, Eckert M, Holman GM, Nachman RJ, Wang Y, et al. Isolation of periviscerokinin-2 from the abdominal perisymphathetic organs of the American cockroach, *Periplaneta americana*. *Peptides* 1998;19:801–9.
- [24] Predel R, Russell WK, Tichy SE, Russell DH, Nachman RJ. Peptidomics of CNS-associated neurohemal systems of adult *Drosophila melanogaster*: a mass spectrometric survey of peptides from individual flies. *Peptides* 2003;24:1487–91.
- [25] Predel R, Wegener C, Russell WK, Tichy SE, Russell DH, Nachman RJ. Peptidomics of CNS-associated neurohemal systems of adult *Drosophila melanogaster*: a mass spectrometric survey of peptides from individual flies. *J Comp Neurol* 2004;474:379–92.
- [26] Riehl MA, Garczynski SF, Crim JW, Hill CA, Brown MR. Neuropeptides and peptide hormones in *Anopheles gambiae*. *Science* 2002;298:172–5.
- [27] Vanden Broeck J. Neuropeptides and their precursors in the fruitfly, *Drosophila melanogaster*. *Peptides* 2001;22:241–54.

- 
- [28] Verleyen P, Huybrechts J, Sas F, Clynen E, Baggerman G, DeLoof A, et al. Neuropeptidomics of the grey flesh fly, *Neobellieria bullata*. *Biochem Biophys Res Commun* 2004;316:763–70.
- [29] Wegener C, Herbert Z, Eckert M, Predel R. The periviscerokinin (PVK) peptide family in insects: evidence for the inclusion of CAP(2b) as a PVK family member. *Peptides* 2002;23:605–11.